NOTA INFORMATIVA

INFORMATIVE NOTE

cod. 393020 mod. Nitrile

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Read this briefing note through carefully prior to use and prior to any maintenance work. The purpose of the information contained herein is to assist and direct the user in choosing and using the equipment. Neither the manufacturer nor the distributor can be held liable for an incorrect use of the equipment. This briefing note must be kept for the duration of the equipment.
Glove model: **cod. 393020 mod. Nitrile**.

Trademark of distributor

The CE 0465 mark means that these gloves are personal protective equipment conforming to the essential health and safety requirements set forth in the directive 89/686/EEC and have been certified by the notified body: **A.N.C.I. Servizi Srl - Sezione C.I.M.A.C. - C.so G. Brodolini, 19 - 27029 Vigevano (PV), Italia, n° 0465**. In addition, the equipment is subject to the control of CE quality guarantee of the finish product by the notified body: **A.N.C.I. Servizi Srl - Sezione C.I.M.A.C. - C.so G. Brodolini, 19 - 27029 Vigevano (PV), Italia, n° 0465**.

SIZES

OUR SIZE REFERENCES	SIZES CONTEMPLATED EN 420
S	7
M	8
L	9
XL	10

The dexterity level of the fingers is equivalent to 5 out of 5.
The gloves meet the safety characteristics only if the size is right and worn properly.

RANGE OF USE

The glove is designed to protect the user against contact with chemical products and micro-organisms (intended as fungus or bacteria, not viruses).
The glove can be used also as a CLASS I MEDICAL DEVICE, pursuant to the legal provisions of Legislative Decree 46/97 and Directive 93/42/EEC, for use in the private and professional sector and within the structures of the National Health Service, provided the warnings on use and care indicated on the glove packaging are observed.
Gloves certified according to the EN 455 standard.

PICTOGRAMS AND PROTECTION LEVELS

EN 374

EN 374

PROTECTION AGAINST CHEMICALS	PROTECTION AGAINST MICRO-ORGANISMS
---------------------------------	---------------------------------------

CHEMICAL RESISTANCE TABLE according to EN 374 - 3:

Chemical	Permeation index (0 - 6)	Time of Passage (min)
G Diethylamine	2	30
K Sodium hydroxide 40%	3	60
L Sulphuric acid 96%	2	30

Time of passage [EN 374 - 3]. The time a chemical takes to pass through the glove material. Do not use the gloves for handling the above-mentioned chemical products for longer than the times stated above.

The gloves' performance levels during the following mechanical tests [EN 388] were as follows: abrasion strength: 0;
resistance to cutting by a blade: 0;
ripping strength: 0;
perforation strength: 0

This device can be used as a **CLASS I, NON-STERILE SINGLE-USE MEDICAL DEVICE** (Rule 5, first section of Annex IX of Directive 93/42/EEC) for medical use, from exploration and hygiene, to generic use.
CND T01020299

COMPLIANT WITH: REGULATION (EU) N. 10/2011 - REGULATION (EC) N. 1935/2004 DIRECTIVES 85/572/EEC AND SUBSEQUENT AMENDMENTS.

Ethanol 20%	Acetic acid 3%	Vegetable oil	Polyphenyl oxide (MPP0)
Yes	No	No	No
1h/40°C	-	-	-

Since the performance levels are based on lab test results, they may differ from the actual conitions at the workplace. Consequently it is advisable to conduct a preliminary test of the gloves to see if they are suitable for their intended use.

The protective characteristics refer to new equipment, in good condition and that has never been subject to cleaning treatments.

WARNINGS

Prior to use, make sure the gloves are in good condition by which we mean no cuts, holes or cracks, etc. If the gloves are not in perfect condition, replace them immediately. The gloves must be used solely for the risks contemplated in this briefing note.
If used as a NON-STERILE medical device, its use must be transient if used invasively, in relation to the natural orifices of the body. It MUST NOT be used in contact with open wounds, or for surgery or inside unnatural orifices of the body or wherever the use of a STERILE device is requested.

STORING

The gloves are packed in cardboard boxes.
Keep the gloves in their original packaging in a clean, dry place, away from sources of heat and direct sunlight. The user must always inspect the gloves to ensure they are sound prior to using them.
The love must be used, provided the packaging is perfectly intact and kept in the preservation conditions stated on the packaging itself, **within the indicated expiry date and no later**. Since these are disposable gloves to be used just the once, they must be replaced at the end of their use – they cannot be used a second time.

NOTES

Used gloves can be contaminated with chemicals or biological substances: dispose of them in compliance with the local regulations in force on this subject (dump, incinerator).
The presence of any allergenic substances is not known to the manufacturer at this time: please report any cases found of hypersensitivity or allergic reactions.
In the case of manufacturing defects this equipment will be replaced.
*** In the event of discrepancies between the different translations, the Italian version will exclusively apply.**

Produced by

for:

Registered Community Trade Mark, no 000596585
c/o UAMI - Alicante - Spain
www.boxerline.com - info@boxerline.com

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